

# PepMute™ Plus siRNA Transfection Reagent

----- A Standard Protocol for siRNA  
Transfection of Mammalian Cells

- ☐ 100 µl  
☐ 500 µl  
☐ 1000 µl



5260 Westview Dr, Suite 100  
Frederick MD 21703  
FAX. 301-560-4919  
TEL. 301-330-5966  
Toll Free. 1-(866)-918-6812  
Email: [info@signagen.com](mailto:info@signagen.com)  
Web: [www.signagen.com](http://www.signagen.com)

This product is for laboratory research ONLY and not for diagnostic use

## Introduction:

PepMute™ Plus Reagent is an upgraded version from PepMute™ siRNA transfection reagent. With addition of several pre-screened hydrophobic groups to its peptide backbone, PepMute™ Plus Reagent gains self-assembly capacity when binding nucleic acids, making PepMute™ Plus Reagent a versatile and most powerful gene delivery tool. PepMute™ Plus Reagent have been validated to effectively and reproducibly transfect single siRNA, And co-transfect DNA/siRNA to variety of mammalian cells.

## Important Guidelines for Transfection:

- PepMute™ Plus reagent was formulated as a powerful siRNA delivery tool. For most adherent cell lines and primary cells, siRNA at ~5.0 nM is basically sufficient to obtain up to 90% gene silencing. For hard-to-transfect cells, we recommend using a final siRNA concentration of 50 nM.
- While the standard protocols for siRNA transfection and siRNA/DNA co-transfection are being given below, optimization is often needed for maximal gene silencing.

## PART I. Standard siRNA Transfection of Adherent Cells

### Step I. Preparation of Working Solution of PepMute™ Transfection Buffer:

PepMute™ Transfection Buffer (5x ) is provided as 5x concentrated stock solution. To make working solution, dilute one part of the stock solution with 4 parts of ddH<sub>2</sub>O into a sterile bottle. The working solution is stable at 4 °C~RT for 12 months.

### Step II. Cell Seeding:

Cells should be plated 18 to 24 hours prior to transfection so that the monolayer cell density reaches to the optimal ~50% confluency at the time of transfection. Complete culture medium with serum and antibiotics is freshly added to each well 30~60 minutes before transfection.

**Note:** PepMute™ Plus reagent is NOT interfered by serum and antibiotics, therefore serum and antibiotic containing medium can be used during the entire experiment.

**Table 1. A Guideline for siRNA transfection per cell culture vessel**

| Culture Dish     | Growth Medium (ml) | Transfection Buffer (µL) | siRNA (pmoles) Final 5.0 or 50 nM | PepMute™ Plus (µL) |
|------------------|--------------------|--------------------------|-----------------------------------|--------------------|
| 24-well          | 0.5                | 50                       | 2.5 / <b>25</b>                   | 1.2 ~ <b>2.0</b>   |
| 12-well          | 0.75               | 75                       | 3.75 / <b>38</b>                  | 2.0 ~ <b>3.3</b>   |
| 6-well           | 1.0                | 100                      | 5.0 / <b>50</b>                   | 2.4 ~ <b>4.0</b>   |
| 60 mm            | 3.0                | 300                      | 15 / <b>150</b>                   | 7.2 ~ <b>12</b>    |
| 10 cm / Flask 75 | 8.0                | 800                      | 40 / <b>400</b>                   | 20 ~ <b>33</b>     |

## Step III. siRNA Transfection Protocol:

For optimal siRNA-mediated silencing, we recommend using 1~100 nM siRNA. As a starting point, we recommend using 5.0 nM siRNA which usually gives satisfactory silencing result for most adherent cell lines or primary cells. For hard-to-transfection cells, we recommend using a final siRNA concentration of 50 nM. (bold & underlined in **Table 1**).

The following conditions are given per well in a 6 well plate. For other culture format, please refer to **Table 1**.

- For each well, add 1.0 ml of complete medium with serum and antibiotics freshly 30~60 minutes before transfection.
- Dilute 5.0 or **50** pmoles siRNA (final concentration of 5.0 or **50** nM respectively per well) into 100 µl of working solution of PepMute™ Transfection Buffer prepared in **Step I**. Pipette up and down to mix.

**Note: For maximum gene silencing, dilute siRNA and PepMute™ Plus reagent with PepMute™ Transfection Buffer (1x).**

**We strongly suggest reconstituting siRNA stock solution at 10.0 µM, so add 0.5 or 5.0 µl siRNA stock solution per well of 6-well plate to make final 5.0 and 50 nM siRNA respectively.**

- Add 2.4 µl or 4.0 µl (for hard-to-transfect cells, bold and underlined in **Table 1**) PepMute™ Plus reagent, mix by pipetting up and down.
- Incubate for ~15 minutes at RT to let transfection complex form.
- Note: Never keep the complex longer than 30 minutes.**
- Add the transfection mix to the cells drop wise. Gently rock the plate back and forth and return the plate to CO<sub>2</sub> incubator.
- Gene silencing is usually measured 24~78 hours post transfection.

## PART II. A Standard Protocol for DNA/siRNA Co-transfection

### Step I. Preparation of Working Solution of PepMute™ Transfection Buffer:

PepMute™ Transfection Buffer (5x ) is provided as 5x concentrated stock solution. To make working solution, dilute one part of the stock solution with 4 parts of ddH<sub>2</sub>O into a sterile bottle. The working solution is stable at 4 °C~RT for 12 months.

# PepMute™ Plus siRNA Transfection Reagent

----- A Standard Protocol for siRNA/DNA  
Co-transfection of Mammalian Cells

- ☐ 100 µl
- ☐ 500 µl
- ☐ 1000 µl



5260 Westview Dr, Suite 100  
Frederick MD 21703  
FAX. 301-560-4919  
TEL. 301-330-5966  
Toll Free. 1-(866)-918-6812  
Email: [info@signagen.com](mailto:info@signagen.com) Web:  
[www.signagen.com](http://www.signagen.com)

This product is for laboratory research ONLY and not for diagnostic use

## Step II. Cell Seeding:

Cells should be plated 18 to 24 hours prior to transfection so that the monolayer cell density reaches to the optimal ~70% confluency at the time of transfection. Complete culture medium with serum and antibiotics is freshly added to each well 30~60 min before transfection.

**Note:** PepMute™ Plus reagent is NOT interfered by serum and antibiotics, therefore serum and antibiotic containing medium can be used during the entire experiment.

**Table 2. A Guideline for DNA & siRNA Co-transfection Per Cell Culture Vessel**

| Culture Dish     | Growth Medium (ml) | Transfection Buffer (µL) | Plasmid DNA (µg) | siRNA (pmoles) Final 5.0 nM | PepMute™ Plus (µL) |
|------------------|--------------------|--------------------------|------------------|-----------------------------|--------------------|
| 24-well          | 0.5                | 50                       | 0.25             | 2.5                         | 1.5                |
| 12-well          | 0.75               | 75                       | 0.375            | 3.25                        | 2.25               |
| 6-well           | 1.0                | 100                      | 0.5              | 5.0                         | 3                  |
| 60 mm            | 3.0                | 300                      | 1.5              | 15                          | 9                  |
| 10 cm / flask 75 | 8.0                | 800                      | 4.0              | 40                          | 24                 |

## Step III. DNA & siRNA co-transfection protocol:

For DNA/siRNA co-transfection experiment, we recommend using 0.25~0.5 µg DNA and 1~20 nM siRNA per well in a 6-well plate.

As a starting point, we recommend using 0.5 µg DNA and 5.0 pmoles siRNA (final concentration 5.0 nM) per well of a 6-well plate which usually give satisfactory knockdown effect.

The following conditions are given per well of a 6 well plate. For other culture format, please refer to **Table 2**.

- For each well, add 1.0 ml of complete medium with serum and antibiotics freshly 30~60 minutes before transfection.
- Dilute 0.5 µg DNA and 5.0 pmoles siRNA (final 5.0 nM) into 100 µl of working solution of PepMute™ Transfection Buffer. Vortex to mix followed by brief spin to bring drops to the bottom of the tube.

**Note:** For optimal transfection efficiency and maximum gene silencing, PepMute™ Transfection Buffer is a must for diluting siRNA/DNA and PepMute™ reagent. We strongly suggest preparing siRNA stock solution at 5.0 µM, so add 1.0 µl siRNA stock solution per well of 6-well plate to make final 5.0 nM of siRNA.

- Add 3 µl PepMute™ Plus reagent immediately, mix by pipetting up and down.
- Incubate for ~15 min at RT to let transfection complex form.
- Note:** Never keep the complex longer than 30 minutes.
- Add the transfection complex to the cells drop wise.

- Gently rock the plate back and forth and return the plate to the incubator.
- Replace transfection medium by cell growth medium ~5 hours after transfection when necessary.
- Gene silencing is usually measured 24~72 hours post transfection.

**Storage:** PepMute™ Plus siRNA Transfection Reagent is stable for up to 12 months at 4 °C. This item shipped at ambient temperature